

Tel: 1-888-762-2568 Fax: 1-510-783-5386 Em ail: info@biochain.com

Certificate of Analysis

Product Name: Membrane Protein Lysate - Hun	nan Tum	or Tissue: Sple	een		
Catalog No.: P3235246	Lot No.: B103045				
Species: ■ Human □ Mouse □ Rat □ Monl	key (Rh)	□ Guinea Pig	□ Porcine		
□ Bovine □ Hamster □ Dog □ Monk	ey (Cy)	□ Rabbit	□ Plant		
Tissue Type: ☐ Normal ☐ Adult ☐ Fetal ■	Tumor	□ Disease	☐ Cell Line		
Donor Information:					
Male: <u>40</u> year(s) old Pool o	f <u> </u>	donor(s)			
Female: year(s) old Pool o	f	_donor(s)			
Pathological Diagnosis: Cappillary Hemangior	na				
Packaging: 0.1 mg per vial; Concentration: 2 n					
Quality Control: Protein lysate is verified on SDS-	PAGE co	omparing to star	ndard protein lysate.		
MW (kDa)					
250					
130					
95	10.7				
72 55					
36 28					
17					

Protein lysate can be validated by Western Blot using specific biomarker: GAPDH for total protein, cytoplasmic and universal total protein; Histone for nuclear protein; Na/K ATPase for membrane protein; Vimentin for cytoskeleton protein.

Use s: Western Blot (10 μg to 20 μg per lane is recommended for a mini gel), Immunoprecipitation, Electrophoresis, Enzymatic activity analysis, protein-protein interaction, gel shift assay, tissue specific expression, and mass spec. analysis.

	u	f	f	Δ	r	•
ם	u	ı	ı	C	ı	•

<u>Location</u>	Buffer Components*
Total protein	HEPES (pH7.9), KCl, Sucrose, Sodium deoxycholate, and NP-40
Membrane protein	The same as above.
Nuclear protein	HEPES (pH7.9), NaCl
Cytoplasmic protein	HEPES (pH7.9), KCI, and Sucrose.
Cytoskeleton protein	PIPES (pH 6.8), NaCl, SDS, and Sucrose

Note: *Each buffer contains, MagCl₂, EDTA, Glycerol, and a cocktail of protease inhibitors

Storage: Store at 2-8°C for continuous use within a day. For extended storage, freeze working aliquots at -70°C. Repeated freezing and thawing is not recommended. Under proper storage conditions the shelf life is six months from the date of receipt.

Shipping: Dry Ice

User's Note: Products are intended for laboratory research purposes only and should be used by qualified personnel only.

APPROVED BY: